


DESKTOP PORTAL


- 1 Select a Patient.
- 2 Click Test Results.
- 3 Click result to view in bottom window.
- 4 Click Images link. The Images will open in a new window for review.














The desktop portal interface shows a patient list on the left. A patient is selected, and the 'Test Results' tab is active. A table of test results is displayed, with one result highlighted. A pop-up window titled 'Universal Viewer Zero Footprint - Google Chrome' shows a chest X-ray image. The 'Images' link in the test result details is highlighted.

Date/Time	Test (10 of 28)	Status	Source
03/17/21 06:35	CHEST PA & LAT	Signed	Trident Regional M
03/12/21 18:33	PORTABLE CHEST	Signed	Trident Regional
08/31/20 08:10	CHEST PA & LAT	Signed	Trident Regional
08/13/20 06:23	PORTABLE CHEST	Signed	Trident Regional
08/12/20 19:15	PORTABLE CHEST	Signed	Trident Regional
08/10/20 18:00	PORTABLE CHEST	Signed	Trident Regional
04/15/20 06:30	CHEST PA & LAT	Signed	Trident Regional
04/14/20 06:33	CHEST PA & LAT	Signed	Trident Regional
04/09/18 11:50	CHEST PA & LAT	Signed	Grand Strand Med
10/17/16 10:20	CHEST PA & LAT	Signed	Grand Strand Med

MOBILE

- 1 From the Patient Chart screen, tap Test Results
- 2 Tap the name of the result to view the result or the  icon to open the image.
- 3 Use the toolbar at the bottom to navigate through the image.



-  **Series:** Tap to open & select series to view
-  **Window Leveling:** Tap to enter & exit window level adjust mode
-  **Swipe up & down** to adjust window levels & contrast.
-  **Leveling Preset:** Turn on/off Invert or Recommended.
-  **Info Overlay:** Turn patient info overlay on or off
-  **Cine Mode:** Scroll up & down to cine mode through image slices
-  **Zoom:** Pinch in & out to zoom in & out
-  **Pan:** Use two fingers to move to different parts of image on screen