

DESKTOP PORTAL

1 Select a Patient.

2 Click Clinical Notes.

The screenshot shows the PatientKeeper desktop interface. On the left is a 'Patient List' with columns for patient name, age, gender, and LOS. The patient 'CPOES, TEST3' (46Y Female LOS:76D) is selected. The main area shows 'Clinical Notes' for this patient. A table lists notes with columns for Date/Time, Note Type, Author, and Source. One note is highlighted: '03/12/21 20:38 EKG - CONFIRMED' by 'SCMUSE EKG'. Below this, a window titled 'MUSE(R) Report For' displays an EKG waveform on pink paper with leads I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, and V6. The report includes 'Referred by: SELF REFERRED' and 'Confirmed By:'.

3 EKGs will display as the following Note Types in the Clinical Notes list. Click to view them in the bottom window:

- EKG – CONFIRMED (read)
- EKG – UNCONFIRMED (pending)

Tip: Click the Column Header Names to sort the list based on preference.

MOBILE

1 Tap Patient to select.

2 Tap Clinical Notes.

3 Find & tap EKG – CONFIRMED or EKG – UNCONFIRMED.

The mobile app navigation is shown in four steps:
 

- Step 1:** The 'Z\_Test Patients List' screen shows a list of patients. 'CPOES, TEST3' is highlighted with a red box.
- Step 2:** The 'Clinical Notes' screen for 'CPOES, TEST3' is shown. The 'Clinical Notes' icon is highlighted with a red box.
- Step 3:** The 'Clinical Notes (45)' list for 'CPOES, TEST3' is shown. The 'EKG - UNCONFIRMED' note is highlighted with a red box.
- Step 4:** The EKG waveform is displayed in a window. A hand icon with a double-headed arrow indicates a pinch-to-zoom gesture on the waveform.

4 The EKG will display in the window.

Pinch on the screen to zoom in and out.